
1. INTENDED USE OF SPINESCAN3D

The SpineScan3D is intended for screening the spine by capturing its tilts and degree of rotations along the spine length. To operate the SpineScan3D, a mobile phone with Android or iPhone platforms should be in place with a mobile application (App) "SpineScan3D" being installed to serve as the control interface. User must follow the instruction in this document.

The SpineScan3D is a system that includes the hardware, software, network, App, workflow and algorithm, which are developed for Research Purpose Only and NOT For Sale. User must NOT use the SpineScan3D for purposes other than stipulated and are bound to the Terms and Conditions.

2. DESCRIPTION OF SPINESCAN3D DEVICE

The SpineScan3D device has several functional features as described below:

5. PROCEDURES FOR USE

Installation of the Android App "SpineScan3D"

5.1. The App supports Android version 4.3 or above.

5.2. Install the software via either:

5.2.1. Google Play Store, or;

5.2.2. Download the App from the download section of the company's website (www.spinescan3d.com) using the browser in the Android phone.

5.2.2.1. Make sure that your phone allows installation from "Unknown Sources" in the default setting.

5.2.2.2. When download is complete, tap on the file to start installation.

5.2.2.3. It is recommended to disable the "Unknown Sources" after the installation.

Installation of the iOS App "SpineScan3D"

5.3. The App supports iOS version 8.0 or above.

5.4. Install the software via App Store.

5.4.1. Open the App Store on your iPhone to browse SpineScan3D App.

5.4.2. Tap on the App to view the App's product page.

5.4.3. Tap "Get" to get the App.

5.4.4. Might need to enter the password to complete your download and installation.

6. REGISTRATION BEFORE USE

For first time User, registration is required before use.

6.1. Launch the App by clicking the icon "SpineScan3D" on the phone.

6.2. Press "Register" to start the registration procedure.

6.3. Enter the following information:

6.3.1. Email Address

6.3.2. Password

6.3.3. Country of operation

6.3.4. Role (select "Doctor" or "Nurse" from the pull-down list)

7. ROUTINE OPERATION PROCEDURE

7.1. Ensure 3 x AAA batteries are inserted as shown in the device.

7.2. Disinfect the wheels with 70% alcohol swab/wipe.

7.3. Mount the phone on the extendable holder of the device. The arrow printed on the holder shows the extending direction. Pull the holder gently according to the direction displayed.

Caution: Avoid the side buttons of the phone being pressed by the holder. Adjust the phone position when it occurs.

7.4. Turn on the device by pressing the "Power" button.

7.5. The Red-Green LED flickers in a normal startup. If battery is low, the Red LED flickers and battery replacement is needed.

7.6. Launch the "SpineScan3D" App on the phone.

7.7. Login to the registered account.

7.8. Ensure the Bluetooth function of the phone is on. Click "Open" and then "Allow" if the Bluetooth function is off.

7.9. Click "Connect" to connect the device with the phone via Bluetooth.

7.10. Select the device to pair from the device list. The device name has prefix of "SpineS..." and is labelled on the device body.

7.11. After the connection is established, Press "Scan" on the phone to start the scanning procedure.

7.12. Input Subject data:

7.12.1. For a first-time Subject;

7.12.1.1. Press "+ New Subject" at the bottom on the "Select Subject" screen.

7.12.1.2. Input the Subject related information as requested.

7.12.1.3. Press "Save" to create the Subject record.

7.12.2. For a re-visiting Subject;

7.12.2.1. Search by entering the ID on the search field.

7.12.2.2. Select the Subject from the search result list.

7.13. User can opt to take a photo record of the Subject's back for the scan. Local regulatory requirements for electronic records storage apply.

7.14. Press "Start Scanning" to proceed.

7.15. The Subject is required to stand in the **Forward-Bending Posture** as shown in the App interface on the phone.

7.16. Place the device at the buttock cleft of the Subject. Ensure the Scanning direction is correct.

7.17. Press the "SCAN" button on the device for starting the measurement. The LED turns on and stays on GREEN.

3. PRECAUTIONS

● User is required to read all the instruction and watch the instruction video in the App before use. The instruction video can be viewed anytime from the App "SpineScan3D".

● Do not immerse the device into any liquid.

● Operating temperature and humidity is in the range of 18-35 °C and 20-80 %, respectively.

● Do not press too hard against the spine.

● Stop applying the device on the Subject when the Subject feels uncomfortable or pain.

● Do not use the device when it is damaged or defective.

● Do not use the device if the two pairs of wheels are not in alignment or are not leveled if laid resting on horizontal plane.

● User should be cautions if the Subject:

- has spine fracture or a history of spine fracture;

- has wounds on the skin of the application area.

● Only use the device in environments where wireless Bluetooth transmission is permitted.

● For Subject carrying implanted medical devices, always consult their healthcare providers before use.

4. SPECIFICATIONS

Measurement range (Tilt Angle)	90 to 90 deg
Measurement resolution (Tilt Angle)	0.01 deg
Measurement resolution (Change of angular speed)	0.5 deg/s
Displacement measurement resolution	0.5 cm
Connectivity	Bluetooth v4.0 BLE
Power (battery not included)	Alkaline 1.5 V (AAA Size) x 3
Dimension	200 mm (W) x 70 mm (H) x 82 mm (D)
Weight	325 g (excluding battery)
Supported phone width	65-78 mm

6.4. Read and click the box for accepting the Terms and Conditions.

6.5. Press "Register" to complete the registration.

6.6. An email notification with a link to activate the account will be sent to the registered email address.

6.7. Click the link in the email to activate the account.

6.8. Registered User can login to the App homepage after successful registration.

7.18. Starting from the level of the iliac crest, move the device steadily along the spinous process of the vertebral column of the Subject until it reaches C7 of the neck. Ensure the encoder wheels have good contact on the Subject's skin. It is recommended that the spinous process should be marked to ensure the scanning is along the spinous process.

7.19. Real-time tilt reading and the incremental spine distance are displayed. Warning message(s) prompt/prompts if unusual measurement(s) is/are detected.

7.20. Press the "SCAN" button again on the device for completing the measurement. The GREEN LED flickers.

7.21. If unusual measurement(s) is/are detected during the scan, warning message(s) is/are prompted. The User can choose "Retry" to repeat the scan from step 7.15 or "Continue" to save the record and proceed to the next step.

7.22. The App shows the **Standing Posture** for Subject to follow.

7.23. Place the device at the neck of the Subject. Ensure the Scanning direction is correct and prevent initial tilt.

7.24. Press the "SCAN" button on the device to start the measurement. The LED turns on and stays on GREEN.

7.25. Move down the device steadily along the vertebral column of the Subject until it reaches the buttock cleft. Ensure the encoder wheels have good contact on the Subject's skin.

7.26. Press the "SCAN" button again on the device for completing the measurement. The GREEN LED flickers.

7.27. If unusual measurement(s) is/are detected during the scan, warning message(s) is/are prompted. The User can choose "Retry" to repeat the scan from step 7.22 or "Continue" to save the record and complete the measurement.

7.28. A bending scan summary, ("results from **Forward-Bending Posture**" only), is displayed. Press "Back" to go back to home screen of the App.

7.29. After scanning, Subject can opt to fill in a questionnaire regarding the health of his/her spine.

7.30. Disinfect the wheels of the device with 70% alcohol swab/wipe before performing the scan on the next Subject.

7.31. After use, turn off the device, detach the phone and store the device in a cool and dry place.

Data Synchronization

The data will be synchronized automatically when the mobile phone is connected to the Internet with the "SpineScan3D" mobile application launched. Ensure that the mobile phone can be connected to the Internet 1) during or 2) no longer than 1 week after the measurement. This is critical for keeping the data in cloud database updated.

● Warning message will pop up if data has not been synchronized for more than one week.

8. INDICATIONS OF DEVICE STATUS


LED Indicator	Information
Flickering Red	Battery is low in voltage
Flickering Red-Green	Fail in establishing Bluetooth connection
Flickering Green	Succeed in establishing Bluetooth connection
Constant Green	The scan button is pressed and the measurement is in progress

Warning messages will pop up on the interface of the phone through

the "SpineScan3D" App (Refer to Section 10.2.2).

10. WARNING AND LABELS

10.1. SpineScan3D device:

11. REPAIR AND MAINTENANCE

User must **NOT** by his/her own self or incur any other third party to repair or modify the device. If the device or the mobile App is not operating properly, please contact Avalon SpineCare (HK) Limited.

12. DISPOSAL

Please contact the Avalon SpineCare (HK) Limited for procedures to return the device.

13. TERMS AND CONDITIONS

User agrees that the device is licensed to the User for the use as stipulated in this Instruction For Use only. Furthermore, User acknowledges this device is for Research Purpose Only and NOT For Sale, and the use of the device does not amount to any medical treatment or diagnosis or professional advices of any disease or abnormality of the Subjects. User agrees that he or she will use the device in accordance to the local regulation applicable to the device.

All the proprietary right of the device, including but not limited to any technology, patent, copyright, trademarks, designs and trade secrets, remains in Avalon SpineCare (HK) Limited. In any circumstances, User must NOT, either by his or her own self or with any third parties, sell, assign, license, transfer, distribute, modify, rent, lease, loan, replicate, manufacture, reverse engineer the device, including the hardware, the software and related documents, in whole or in any part, without written permission from Avalon SpineCare (HK) Limited.

User expressly acknowledges and agrees that use of the device is at User's sole risk. The device, including the hardware, the software and the related documents are provided "as is" and without warranty of any kind. Avalon SpineCare (HK) Limited expressly disclaims all warranties, expressed or implied, including but not limited to the fitness of a particular purpose. Avalon SpineCare (HK) Limited does not warrant that the functions contained in the device, including the hardware and software, will be uninterrupted or error-free, or that defects will be corrected. Furthermore, Avalon SpineCare (HK) Limited does not warrant or make any representations regarding the results of the use of the devices in terms of the correctness, accuracy, reliability, or otherwise. No oral or written information or advice given by Avalon SpineCare (HK) Limited, including but not limited to the directors, employees nor representatives, create a warranty or in any way increase the scope of this warranty. Should the device prove defective, user contacts Avalon SpineCare (HK) Limited for arrangements in necessary servicing, repair or correction. Avalon SpineCare (HK) Limited does not warrant or guarantee that the features, benefits and functionality of the device will remain the same in future updates. User understands and agrees that any feature may change at the sole discretion of Avalon SpineCare (HK) Limited and that the use of any feature or function of the device may be discontinued or modified at the sole discretion of Avalon SpineCare (HK) Limited without notice to User.

Avalon SpineCare (HK) Limited shall not be held liable for any damages whatsoever resulting from the use or installation of this product regardless of the cause, including but not limited to defect in material or workmanship, misuse or improper installation of the device, accident regardless of cause, or negligent or faulty design. User hereby waives all liability claims against Avalon SpineCare (HK) Limited, its employees and directors, its partners and affiliates and its distributors. Avalon SpineCare (HK) Limited reserves the right to amend this Instruction For Use from time to time, without notice to the User.

9. STORAGE CONDITIONS

● Store the device in cool and dry place.

● Remove the battery when the device is not used for a long time.

10.2. SpineScan3D mobile application:

10.2.1. In the scan summary, the tilt angle is coloured according to the range of angle.

Default Range [Unsigned Degree]	Colour
0.0-2.9	White
3.0-4.9	Yellow
5.0-14.9	Red
> 15.0	Red

10.2.2. When unusual measurements are detected, a warning message is appended next to the icon in the App interface:

Warning Description	Suggested Action
Initial bending angle out of range (Forward-Bending Posture)	Check Subject's posture and placing position of the device. Retry if necessary.
Forward tilt angle exceed limit	Check Subject's posture and avoid accidental movement during operation. Retry if necessary.
Initial bending angle out of range (Standing Posture)	Check Subject's posture and placing position of the device. Retry if necessary.
Device moving too fast - exceed maximum speed	Retry and slow down the measurement speed.
Device moving too slow - fail to reach minimum speed	Retry and increase measurement speed to ensure continuous measurement.
Measured tilt angle range 3.0-4.9 [unsigned degree]	Check on the Subject's posture and spine. Retry if necessary.
Measured tilt angle range 5.0-14.9 [unsigned degree]	Check on the Subject's posture and spine. Retry if necessary.
Measured tilt angle > 15.0 unsigned degree	Check on the Subject's posture and spine. Retry if necessary.

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